



# Rabbit anti-Human GLB1L2 Polyclonal antibody (CABT-B9627)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	Recombinant Protein, antigen sequence: RGRVNYGENIDDQRKGLIGNLYLNDSPKLNFRISLDMKKSFFQRFGLDKWSSL
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Affinity purified
Conjugate	Unconjugated
Applications	IHC
Format	Liquid
Concentration	0.1 mg/ml
Size	100 µl
Buffer	40% glycerol and PBS (pH 7.2). 0.02% sodium azide is added as preservative.
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C.

## BACKGROUND

**Keywords** GLB1L2; galactosidase, beta 1-like 2; MST114; MSTP114; beta-galactosidase-1-like protein 2;

# GENE INFORMATION

Entrez Gene ID [89944](#)

---

UniProt ID [A0A024R3N9](#)

---