



Rabbit anti-Human ST6GALNAC3 Polyclonal antibody (CABT-B9629)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant Protein, antigen sequence: EVNFPLLLNCFGQPGTKWIPFSYTYRRPLRTHYGYINVKTQEPLQLDCDL
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Affinity purified
Conjugate	Unconjugated
Applications	IHC, WB, ICC-IF
Format	Liquid
Concentration	0.05 mg/ml
Size	100 µl
Buffer	40% glycerol and PBS (pH 7.2). 0.02% sodium azide is added as preservative.
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C.

BACKGROUND

Introduction	ST6GALNAC3 belongs to a family of sialyltransferases that transfer sialic acids from CMP-sialic acid to terminal positions of carbohydrate groups in glycoproteins and glycolipids.
---------------------	---

Keywords ST6GALNAC3; ST6 (alpha-N-acetyl-neuraminy-2,3-beta-galactosyl-1,3)-N-acetylglactosaminide alpha-2,6-sialyltransferase 3; STY; SIAT7C; PRO7177; ST6GALNACIII; alpha-N-acetylglactosaminide alpha-2,6-sialyltransferase 3; SIAT7-C; ST6GALNAC III; sialyltransferase 7C; GalNAc alpha-2,6-sialyltransferase III; alpha-N-acetylglactosaminide alpha-2,6-sialyltransferase III; ST6 (alpha-N-acetyl-neuraminy-2,3-beta-galactosyl-1,3)-N-acetylglactosaminide alpha-2,6-sialyltran; sialyltransferase 7 ((alpha-N-acetylneuraminy-2,3-beta-galactosyl-1,3)-N-acetyl galactosaminide alpha-2,6-sialyltransferase) C;

GENE INFORMATION

Entrez Gene ID [256435](#)

UniProt ID [Q8NDV1](#)
