



# Rabbit anti-Human ARPP19 Polyclonal antibody (CABT-B9631)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Recombinant Protein, antigen sequence: MSAEVPEAASAEQKEMEDKVTSPKAEAAK
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Affinity purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	IHC
<b>Format</b>	Liquid
<b>Concentration</b>	0.1 mg/ml
<b>Size</b>	100 µl
<b>Buffer</b>	40% glycerol and PBS (pH 7.2). 0.02% sodium azide is added as preservative.
<b>Storage</b>	Store at +4°C for short term storage. Long time storage is recommended at -20°C.

## BACKGROUND

<b>Introduction</b>	The 19-kD cAMP-regulated phosphoprotein plays a role in regulating mitosis by inhibiting protein phosphatase-2A.
<b>Keywords</b>	ARPP19; cAMP-regulated phosphoprotein, 19kDa; ENSAL; ARPP16; ARPP-16; ARPP-19;

## GENE INFORMATION

Entrez Gene ID [10776](#)

---

UniProt ID [A0A024R5U3](#)

---