



Rabbit anti-Human ABHD16B Polyclonal antibody (CABT-B9634)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant Protein, antigen sequence: VVMAREGRAVVTRWLRAGSLAQEAIFYARYRVDEDWCLALLRSYRARCEEELEGEEALGP HGPAFPWLVGQ
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Affinity purified
Conjugate	Unconjugated
Applications	IHC
Format	Liquid
Concentration	0.1 mg/ml
Size	100 µl
Buffer	40% glycerol and PBS (pH 7.2). 0.02% sodium azide is added as preservative.
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C.

BACKGROUND

Keywords	ABHD16B; abhydrolase domain containing 16B; C20orf135; dJ591C20.1; abhydrolase domain-containing protein 16B;
-----------------	---

GENE INFORMATION

Entrez Gene ID [140701](#)

UniProt ID [Q9H3Z7](#)
