



Rabbit anti-Human GOLST1B Polyclonal antibody (CABT-B9640)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant Protein, antigen sequence: RRVPVLGSLNLPGIRSFVDKVGESN
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Affinity purified
Conjugate	Unconjugated
Applications	IHC
Format	Liquid
Concentration	0.1 mg/ml
Size	100 µl
Buffer	40% glycerol and PBS (pH 7.2). 0.02% sodium azide is added as preservative.
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C.

BACKGROUND

Keywords GOLST1B; golgi transport 1B; GCT2; GOT1; GOT1B; CGI-141; YMR292W; vesicle transport protein GOT1B; hGOT1a; germ cell tumor 2; putative NF-kappa-B-activating protein 470;

GENE INFORMATION

Entrez Gene ID [51026](#)

UniProt ID [Q9Y3E0](#)
