



Rabbit anti-Human C2orf27B polyclonal antibody (CABT-B9644)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant Protein, antigen sequence: FPVPAVEVPVPAPGADSPPGTALEEEAPEPSCRCPGTAQDQPSEELPDFMAPPVEPPASA LELKVWLELEVAERGGQHSSSQQLPHCSQSWAQWKLWRQRPGFIAIWAPLPHWRGTSLIQQ SSSPAAEGPAATAAGAVCLPAGGAGEQEKEPVSRGSSRSSCSQRRPPPPGMEVCPQLGIW AICP(26-209 aa encoded by BC043584)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	WB, ELISA
Positive Control	Transfected HEK-293 cells
Format	Liquid
Size	100 µl
Buffer	PBS with 0.1% sodium azide and 50% glycerol pH 7.3.
Storage	Store at -20°C.

BACKGROUND

Keywords

C2ORF27B; chromosome 2 open reading frame 27B; uncharacterized protein LOC408029;

GENE INFORMATION

Entrez Gene ID[408029](#)**UniProt ID**[Q580R0](#)
