



Rabbit anti-Human C1QTNF6 polyclonal antibody (CABT-B9646)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Isotype	IgG
Source/Host	Rabbit
Species Reactivity	human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	WB, ELISA
Molecular Weight	55 kDa
Positive Control	HeLa cells
Format	Liquid
Size	100 µl
Buffer	PBS with 0.1% sodium azide and 50% glycerol pH 7.3.
Storage	Store at -20°C.

BACKGROUND

Keywords	C1QTNF6; C1q and tumor necrosis factor related protein 6; CTFP6; CTRP6; ZACRP6; complement C1q tumor necrosis factor-related protein 6; complement-c1q tumor necrosis factor-related protein 6;
----------	---

GENE INFORMATION

Entrez Gene ID [114904](#)

UniProt ID [A0A024R1J0](#)
