



Rabbit anti-Human DUS3L polyclonal antibody (CABT-B9650)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant Protein, antigen sequence:MAEGTAEAPLENGGGGDSGAGALERGVAPIKRQYLTTKEQFHQFLEAKGQE KTCRETELDIRGKLYLAPLTTCGNLPFRRICKRFGADVTCGEMAVCTNLLQGQMSEWALL KRHQCEDIFGVQLEGAFPDTMTKCAELLSRTVEVDFVDINVGCPIDLVIKKGGGCALMNR STKFQQIVRGMNQVLDVPLTVKIRTGVQERVNLAHRLPELRDWGVALVTLHGRSREQRY TKLADWQYIEECVQAASPMPLFGNGDILSFEDANRAMQTGVTGIMIARGALLKPWLFTEI KEQRHWDISSERLDILRDFTNYGLEHWGSDTQGVKTRRFLLEWLSFLCRY(1-343 aa encoded by BC004549)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	human, Mouse, Rat
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	WB, IHC, ELISA
Molecular Weight	72 kDa
Positive Control	MCF7 cells, HeLa cells, Jurkat cells
Format	Liquid
Size	100 µl
Buffer	PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C.

BACKGROUND

Keywords

DUS3L; dihydrouridine synthase 3-like (*S. cerevisiae*); DUS3; tRNA-dihyrouridine(47) synthase [NAD(P)(+)]-like;

GENE INFORMATION

Entrez Gene ID

[56931](#)

UniProt ID

[Q96G46](#)
