



Rabbit anti-Human IFRG15 polyclonal antibody (CABT-B9660)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant Protein, antigen sequence:MFSDNSHCPDCGQQWFPSLELGHWLYQTELVENECYQVFLDRINRADYCPE CYPDNPANRSLVLPWSFPLEWAPQNLTRWTFEKACHPFLGPPPLVRKRIHDSRVAGFNPA LQLILTRTDKTLNKKLGQNK (1-131 aa encoded by BC098348)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	IHC, ELISA
Format	Liquid
Size	100 µl
Buffer	PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Storage	Store at -20°C.

BACKGROUND

Keywords IFRG15; interferon responsive gene 15;

GENE INFORMATION

Entrez Gene ID

[64163](#)
