



Rabbit anti-Human KIAA0895 polyclonal antibody (CABT-B9663)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant Protein, antigen sequence: MVATGWARRSGLHPAPRPLRGAQEALAKNLWREERNRPLILSSLIIKKLHW PEQELAKKSILNAEDSLIIDNKRSISHLSSGVLKDIFTTGTSSYNVLLQSKEEKYHSQK QSSSTYSKRCKRPSKSPNTSRSKDPRRMKALVPVTSSGTWYCLERRPAVFVTSSVSSPVK FTHDISVTGNIGVLPKPKSKVKWCHFSTLPKPKQLSRSEKGGDFSGKKFCILTAIKP TNLEKEKLRFFKSDYTYPQFEYANPALPSVLAKHSHASDRFLKQVKCHFNSGSSFLCTK TIICDLISKNKDI (1-304aa encoded by BC110422)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	human, Mouse
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	WB, ELISA
Molecular Weight	45, 60 kDa
Positive Control	mouse testis tissue, HEK-293 cells
Format	Liquid
Size	100 µl
Buffer	PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Storage	Store at -20°C.

BACKGROUND

Keywords KIAA0895; uncharacterized protein KIAA0895;

GENE INFORMATION

Entrez Gene ID [23366](#)

UniProt ID [Q8NCT3](#)
