



Rabbit anti-Human LRRC20 polyclonal antibody (CABT-B9667)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant Protein, antigen sequence:MLKKMGEAVARVARKVNETVESGSDTLDLAECKLVSFPIGIYKVLNRNVSGQ IHLITLANNELKSLTSKFMTTFSQLRELHLEGNFLHRLPSEVSALQHLKAIDLSRNQFQD FPEQLTALPALETINLEENEIVDVPVEKLAAMPALRSINLRFNPLNAEVRVIAPPLIKFD MLMSPEGARAPLP (1-184 aa encoded by BC024001)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	human, Mouse
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	WB, IHC, IF, ELISA
Molecular Weight	21-25 kDa
Positive Control	NIH/3T3 cells, Y79 cells
Format	Liquid
Size	100 µl
Buffer	PBS with 0.1% sodium azide and 50% glycerol pH 7.3.
Storage	Store at -20°C.

BACKGROUND

Keywords

LRRC20; leucine rich repeat containing 20; leucine-rich repeat-containing protein 20;

GENE INFORMATION

Entrez Gene ID

[55222](#)

UniProt ID

[A0A024QZM2](#)
