



Rabbit anti-Human NRSN2 polyclonal antibody (CABT-B9673)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant Protein, antigen sequence: MMPSCNRSCSCSRGPSVEDGKWYGVRSYLHLFYEDCAGTALSDDPEGPPVL CPRRPWPSLCWKISLSSGTLTLLLLGVAALTGTGYAVPPKLEGIGEGEFLVLDQRAADYNQA LGTCLAGTALCVAAGVLLAICLFWAMIGWLSQDTKAEPLDPEADSHVEVFGDEPEQQLS PIFRNASGQSWFSPASPFGQSSVQTIQPKRDS (1-204 aa encoded by BC001963)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	WB, IHC, ELISA
Molecular Weight	26 kDa
Positive Control	human brain tissue
Format	Liquid
Size	100 µl
Buffer	PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Storage	Store at -20°C.

BACKGROUND

Keywords NRSN2; neurensin 2; C20orf98; dJ1103G7.6; neurensin-2;

GENE INFORMATION

Entrez Gene ID [80023](#)

UniProt ID [Q9GZP1](#)
