



Rabbit anti-Human PSD2 polyclonal antibody (CABT-B9675)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant Protein, antigen sequence:TDLHGHNIGKKMSCQQFIANLDQLNDGQDFAKDLLKTLYNSIKNEKLEWAI DEDELKSLSELVDDKFGTGKKVTRILDGGNPFLDVPQALSATTYKHGVLTRKTHADMD GKRTPRGRRGWKKFYAVLKGITILYLQKDEYRPDKALSEGDLKNAIRVHHALATRASDYSK KSNVLKLTADWRVFLFQAPSKEEMLSWILRINLVAAIFSAPAFPAAVSSMKKFCRPLLP SCTTRLCQEEQLRSHENKLRQLTAELAEHRCHPVERGIKSKEAEEYRLKEHYLTFEKSRY ETYIHLLAMKIKVGSDDLRIEARLATLEGDDPSLRKTHSSPALSQGHVTGSKTTKDATG PDT (418-771aa encoded by BC038233)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	human, Mouse, Rat
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	WB, ELISA
Molecular Weight	85-95 kDa
Positive Control	mouse brain tissue
Format	Liquid
Size	100 µl
Buffer	PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C.

BACKGROUND

Keywords

PSD2; pleckstrin and Sec7 domain containing 2; EFA6C; PH and SEC7 domain-containing protein 2; 6330404E20Rik; pleckstrin homology and SEC7 domain-containing protein 2;

GENE INFORMATION

Entrez Gene ID

[84249](#)

UniProt ID

[Q9BQI7](#)
