



Rabbit anti-Human PLAC8L1 polyclonal antibody (CABT-B9679)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant Protein, antigen sequence: MNWFGSNFFRCPEDLSLLNIYSPLLSHMSEDEHFISNLRGHVPASAVVKQ PVRGASGRTTITAIQTGGGWSTGLFSVCRDRRICFCGLFCPMCLECDIARHYGECLCWP LLPGSTFALRIGTRERHKKIQTLCEDWLAVHCCWAFSICQVARELKMRTSQVYEICAVPM TKDTLV (1-177 aa encoded by BC130588)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	IHC, ELISA
Format	Liquid
Size	100 µl
Buffer	PBS with 0.1% sodium azide and 50% glycerol pH 7.3.
Storage	Store at -20°C.

BACKGROUND

Keywords	PLAC8L1; PLAC8-like 1; PLAC8-like protein 1;
-----------------	--

GENE INFORMATION

Entrez Gene ID [153770](#)

UniProt ID [A1L4L8](#)
