



Rabbit anti-Human RBM12B polyclonal antibody (CABT-B9681)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant Protein, antigen sequence:TGLTIPDGGVHIIGGEIGEAFIIFATDEDARRAISRSGGFIKDSSVELFLS SKAEMQKTIEMKRTRVGRGRPGSGTSGVDSLNSFIESVKEEASNSGYGSSINQDAGFHT NGTGHGNLRPRKTRPLKAENPYLFLRGLPYLVNEDDVRVFFSGLCVDGVIFLKHHHDGRNN GDAIVKFASCVDASGGLKCHRSFMGSRFIEVMQGSEQQWIEFGGNAVKEGDVLRRFEEHS PPRGINDRHFRKRSHSKSPRRTSRSPPLGFYVHLKNLSLSIDERDLRNFFRGTDLTDEQI RFLYKDENRTRYAFVMFKTLKDYNLTALSLHKTVLQYRPVHI (24-355 aa encoded by BC039260)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	human, Mouse, Rat
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	WB, ELISA
Molecular Weight	118 kDa
Positive Control	A431 cells
Format	Liquid
Size	100 µl
Buffer	PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C.

BACKGROUND

Keywords

RBM12B; RNA binding motif protein 12B; MGC:33837; RNA-binding protein 12B; RNA-binding motif protein 12B;

GENE INFORMATION

Entrez Gene ID

[389677](#)

UniProt ID

[Q8IXT5](#)
