



Rabbit anti-Human SFI1 polyclonal antibody (CABT-B9686)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant Protein, antigen sequence:MAQAHPSSLADLGLAFPDFQSPYKIGDPTEISDPPAAPGQPRGTDPMFLSS NTAHSARKQPRRPHFLLEPAQSQRPPQKPQEHGLGMAQPAAPSLTRPFLAEAPTALVPHSP LPGALSSAPGPKQPPTASTGPELLLLPLSSFMPCGAAAPARVSAQRATPRDKPPVPSSLA SVPDPHLLPGDFSATRAGPGLSTAGSLDLEAELEEIQQLLHYQTTKQNLWSCRRQASS LRRWLELNREEPGPEDQEVEQQVQKELEQVEMQIQLLAEELQAQRQPIGACVARIQALRQ ALC (917-1211aa encoded by BC033613)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	human, Mouse
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	WB, IF, ELISA
Molecular Weight	148 kDa
Positive Control	A549 cells, HepG2 cells, mouse lung tissue
Format	Liquid
Size	100 µl
Buffer	PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Storage	Store at -20°C.

BACKGROUND

Keywords

SFI1; Sfi1 homolog, spindle assembly associated (yeast); PISD; hSfi1p; PPP1R139; protein SFI1 homolog; homolog of yeast Sfi1; protein phosphatase 1, regulatory subunit 139;

GENE INFORMATION

Entrez Gene ID

[9814](#)

UniProt ID

[A8K8P3](#)
