



Rabbit anti-Human TMEM71 polyclonal antibody (CABT-B9689)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant Protein, antigen sequence:MYRISQLMSTPVASSRLEREYAGELSPTCIFPSFTCDSLDGYHSFECGSI DPLTGSHYTCRRSPRLLTNGYYIWTEDSLCDKDG NITLNPSQTSVMYKENLVRIFRKKK RICHSFSSLFNLSTSKSWLHGSIFGDINSSPSEDNWLKGTRRLDTHCNGNETRLL (1- 167aa encoded by BC105783)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	IHC, ELISA
Format	Liquid
Size	100 µl
Buffer	PBS with 0.1% sodium azide and 50% glycerol pH 7.3.
Storage	Store at -20°C.

BACKGROUND

Keywords	TMEM71; transmembrane protein 71;
-----------------	-----------------------------------

GENE INFORMATION

Entrez Gene ID [137835](#)

UniProt ID [Q6P5X7](#)
