



Rabbit anti-Human ZDHHC1 polyclonal antibody (CABT-B9697)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant Protein, antigen sequence:KMRPIQEMEFYMRTFRHMRPEPPGQAGPAAVNAKHSRPASPDPTPGRRDCA GPPVQVEWDRKKPLPWRSPLLLLAMWGPQAPPCLCRKRGRGACIKCERLRPRIRRRGLGP PAAAPARRRIPRTPALCTPLALPAPTTRRRQSPWTRFQWRRRAWAAPLWPPRGAGADSPR WRGRRVRPPFS (304-485 aa encoded by BC021908)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	WB, ELISA
Molecular Weight	48-55 kDa
Positive Control	A549 cells
Format	Liquid
Size	100 µl
Buffer	PBS with 0.1% sodium azide and 50% glycerol pH 7.3.
Storage	Store at -20°C.

BACKGROUND

Keywords ZDHHC1; zinc finger, DHHC-type containing 1; ZNF377; C16orf1; HSU90653; probable palmitoyltransferase ZDHHC1; DHHC-1; zinc finger protein 377; zinc finger, DHHC domain containing 1; DHHC-domain-containing cysteine-rich protein; zinc finger DHHC domain-containing protein 1; DHHC domain-containing cysteine-rich protein 1;

GENE INFORMATION

Entrez Gene ID [29800](#)

UniProt ID [Q8WTX9](#)
