



Rabbit anti-Human VSTM1 polyclonal antibody (CABT-B9698)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant Protein, antigen sequence:CLGYEDEKKNEKPPKPSLHAWPSSVVEAESNVTLKCQAHSQNVTFVLRKVN DSGYKQEQSSAENEAEPFTDLKPKDAGRYFCAYKTTASHEWSESSEHLQLVVTDKHDEL EAPSMKTDTRTIFVAI (14-140aa encoded by BC100942)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	WB, FC, ELISA
Molecular Weight	40 kDa
Positive Control	HL-60 cells
Format	Liquid
Size	100 µl
Buffer	PBS with 0.1% sodium azide and 50% glycerol pH 7.3.
Storage	Store at -20°C.

BACKGROUND

Keywords

VSTM1; V-set and transmembrane domain containing 1; SIRL1; SIRL-1; UNQ3033; V-set and transmembrane domain-containing protein 1; LAIR homolog; OSCAR-like transcript-1; signal inhibitory receptor on leukocytes-1;

GENE INFORMATION

Entrez Gene ID[284415](#)**UniProt ID**[D2DJS6](#)