



Rabbit anti-Human ZNF74 polyclonal antibody (CABT-B9699)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant Protein, antigen sequence:RIHSGEKPFKCNECGKAFSSHAYLIVHRRRIHTGEKPFDCSQCWKAFSCHSS LIVHQRIHTGEKPYKCSECGRAFSQNHCLIKHQKIHSGEKSFKCEKCGEMFNWSSHLTEH QRLHSEGKPLAIQFNKHLLSTYYVPGSLLGAGDAGLRDVPIDALDVAKLLCVVPPRAGR NFSLGSKPRN (464-644 aa encoded by BC013395)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	human, Mouse, Rat
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	WB, IF, ELISA
Molecular Weight	64-72 kDa
Positive Control	Jurkat cells, human brain tissue
Format	Liquid
Size	100 µl
Buffer	PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Storage	Store at -20°C.

BACKGROUND

Keywords

ZNF74; zinc finger protein 74; COS52; hZNF7; ZFP520; ZNF520; zinc finger protein 520;

GENE INFORMATION

Entrez Gene ID

[7625](#)

UniProt ID

[Q16587](#)
