



Mouse anti-Human AIDA monoclonal antibody, clone 3F72D22 (CABT-B9736)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	FLJ12806 (AAH15535, 1 a.a. ~ 307 a.a) full length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	3F72D22
Conjugate	Unconjugated
Applications	WB,IF,sELISA,ELISA
Sequence Similarities	MSEVTRSLLRWGASFRRGADFDSWGQLVEAIDEYQILARHLQKEAQAQHNNSEFTEEQK KTIGKIATCLELRSAAALQSTQSQEELKLEDLKKLEPILKNILTYNKEFPFDVQPVPLRRI LAPGEEENLEFEEDEEEGGAGAGSPDSFPARVPGTLLPRLPSEPGMTLLTIRIEKIGLKD AGQCIDPYITVSVKDLNGIDLTPVQDTPVASRKEDTYVHFNVDIELQKHVEKLTGAAIF FEFKHYKPKKRFTST
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	AIDA; axin interactor, dorsalization associated; C1orf80; axin interactor, dorsalization-associated protein; axin interaction partner and dorsalization antagonist;
-----------------	---

GENE INFORMATION

Entrez Gene ID	64853
-----------------------	-----------------------

UniProt ID	Q96BJ3
-------------------	------------------------

Function	protein domain specific binding
-----------------	---------------------------------
