



Mouse anti-Human AKIRIN2 monoclonal antibody, clone 4E0 (CABT-B9742)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	C6orf166 (AAH00764, 1 a.a. ~ 204 a.a) full length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	Human
Clone	4E0
Conjugate	Unconjugated
Applications	WB,IF,sELISA,ELISA
Sequence Similarities	MACGATLKRTLDFDLLSPASPKRRRCAPLSAPTSAAASPLSAAAATAASFSAASPQK YLRMEPSPFQDVSSRLTTEQILYNIKQEYKRMQKRRHLETSFQQTDPCCTSDAQPHAFLL SGPASPGTSSAASSPLKKEQPLFTLRQVGMICERLLKEREKVVREEYEEILNTKLAEQYD AFVKFTHDQIMRRYGEQPASYVS*
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords

AKIRIN2; akirin 2; FBI1; C6orf166; dJ486L4.2; akirin-2; fourteen-three-three beta interactant 1;

GENE INFORMATION

Entrez Gene ID[55122](#)**UniProt ID**[Q9NW35](#)
