



Mouse anti-Human ANGPTL6 monoclonal antibody, clone 2B22 (CABT-B9757)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	ANGPTL6 (NP_114123, 211 a.a. ~ 300 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	Human
Clone	2B22
Conjugate	Unconjugated
Applications	WB,sELISA,ELISA
Sequence Similarities	GSTSDTSRMLDPAPEPQRDQTQRQQEPMASPMAGHPAVPTKPVGPWQDCAEARQAGHEQ SGVYELRVGRHVSVWCEQQLEGGGWTVI*
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	ANGPTL6; angiopoietin-like 6; AGF; ARP5; angiopoietin-related protein 6; angiopoietin-like protein 6; angiopoietin-related protein 5; angiopoietin-related growth factor;
-----------------	---

GENE INFORMATION

Entrez Gene ID	83854
UniProt ID	Q8NI99
Function	receptor binding