



Mouse anti-Human ANKRA2 monoclonal antibody, clone 2E22 (CABT-B9759)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	ANKRA2 (NP_075526, 1 a.a. ~ 99 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	2E22
Conjugate	Unconjugated
Applications	WB,IF,sELISA,ELISA
Sequence Similarities	MDTSTNLDIGAQLVIEECPSTYSLTGMPDIKIEHPLDPNSEECSAQGVAMGMKFILPNRF DMNVCSRFBVKSLNEEDSKNIQDQVNSDLEVASVLFKAE*
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	ANKRA2; ankyrin repeat, family A (RFXANK-like), 2; ANKRA; ankyrin repeat family A protein 2; RFXANK-like protein 2;
-----------------	---

GENE INFORMATION

Entrez Gene ID	57763
UniProt ID	Q9H9E1
Pathway	Signaling events mediated by HDAC Class II, organism-specific biosystem
Function	histone deacetylase binding; low-density lipoprotein particle binding