



# Mouse anti-Human ANKRD37 monoclonal antibody, clone 7F9 (CABT-B9761)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                              |                                                                                                                                                                            |
|------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Immunogen</b>             | ANKRD37 (NP_859077, 1 a.a. ~ 158 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.                                                     |
| <b>Isotype</b>               | IgG2a                                                                                                                                                                      |
| <b>Source/Host</b>           | Mouse                                                                                                                                                                      |
| <b>Species Reactivity</b>    | Human                                                                                                                                                                      |
| <b>Clone</b>                 | 7F9                                                                                                                                                                        |
| <b>Conjugate</b>             | Unconjugated                                                                                                                                                               |
| <b>Applications</b>          | WB,IF,sELISA,ELISA                                                                                                                                                         |
| <b>Sequence Similarities</b> | MLLLDCNPEVDGLKHLLETGASVNAPDPCKQSPVHLAAGSGLACFLLWQLQTGADLNQQ<br>DVLGEAPLHKA AKVGSLECLSLVASDAQIDL CNKNGQTAEDLAWSCGFDP DCAKFLT TIK<br>CMQTIKASEHPDRNDCVAVLRQKRSLGSEVENTSGKRKC |
| <b>Format</b>                | Liquid                                                                                                                                                                     |
| <b>Size</b>                  | 100 µg                                                                                                                                                                     |
| <b>Buffer</b>                | In 1x PBS, pH 7.2                                                                                                                                                          |
| <b>Storage</b>               | Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.                                                                                                   |

## BACKGROUND

**Keywords** ANKRD37; ankyrin repeat domain 37; Lrp2bp; ankyrin repeat domain-containing protein 37;

hLrp2bp; low density lipoprotein receptor-related protein binding protein; low-density lipoprotein receptor-related protein 2-binding protein;

---

## GENE INFORMATION

Entrez Gene ID [353322](#)

---

UniProt ID [Q7Z713](#)

---