



Mouse anti-Human ANKS1A monoclonal antibody, clone 7C22 (CABT-B9762)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	ANKS1 (NP_056060, 1037 a.a. ~ 1133 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	7C22
Conjugate	Unconjugated
Applications	sELISA,ELISA
Sequence Similarities	HVFSTVDVNLTYEIIITLGQAFEVAYQLALQAQKSRATGASAAEMIETKSSKPVPKPRVG VRKSALEPPDMDQDAQSHASVSWVVDPKPDSKRSLs*
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Ankyrin repeat and SAM domain-containing protein 1A is a protein that in humans is encoded by the ANKS1A gene. Mouse monoclonal antibody raised against a partial recombinant ANKS1.
---------------------	--

Keywords

ANKS1A; ankyrin repeat and sterile alpha motif domain containing 1A; ANKS1; ankyrin repeat and SAM domain-containing protein 1A; Odin; ankyrin repeat and SAM domain containing 1;

GENE INFORMATION

Entrez Gene ID

[23294](#)

UniProt ID

[Q92625](#)

Function

protein binding
