



Mouse anti-Human ANTXR1 monoclonal antibody, clone 300D2440 (CABT-B9765)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to amino acids 481-497 of human ANTXR1.
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	300D2440
Conjugate	Unconjugated
Applications	WB, IHC
Format	Liquid
Buffer	In PBS (0.05% BSA, 0.05% sodium azide)
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords ANTXR1; anthrax toxin receptor 1; Tem8; Antrx1; 2310008J16Rik; 2810405N18Rik; tumor endothelial marker 8;

GENE INFORMATION

Entrez Gene ID	84168
UniProt ID	Q9H6X2
Pathway	Cellular roles of Anthrax toxin, organism-specific biosystem;
Function	actin filament binding; collagen binding; metal ion binding; protein binding; receptor activity; transmembrane signaling receptor activity;