



Mouse anti-Human APOL6 monoclonal antibody, clone 2G8 (CABT-B9777)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	APOL6 (NP_085144, 1 a.a. ~ 77 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	Human
Clone	2G8
Conjugate	Unconjugated
Applications	WB,sELISA,ELISA
Sequence Similarities	MDNQAERESEAGVGLQRDEDDAPLCEDVELQDGDLSPEEKIFLREFPRLKEDLKGNIDKL RALADDIDKTHKKFTK*
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	This gene is a member of the apolipoprotein L gene family. The encoded protein is found in the cytoplasm, where it may affect the movement of lipids or allow the binding of lipids to
---------------------	--

organelles.Mouse monoclonal antibody raised against a partial recombinant APOL6.

Keywords	APOL6; apolipoprotein L 6; 2310076O14Rik; apolipoprotein L6;
-----------------	--

GENE INFORMATION

Entrez Gene ID	80830
-----------------------	-----------------------

UniProt ID	Q9BWW8
-------------------	------------------------

Function	lipid binding; lipid transporter activity
-----------------	---
