



Mouse anti-Human ARL14 monoclonal antibody, clone 3D9 (CABT-B9790)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	ARF7 (AAH34354, 1 a.a. ~ 193 a.a) full length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	3D9
Conjugate	Unconjugated
Applications	sELISA, ELISA
Sequence Similarities	MGSLGSKNPQTKQAQVLLGLDSAGKSTLLYKLKLAKDITTIPTIGFNVEMIELERNFSL TVWDVGGQEKMRVTWGCYCENTDGLVYVVDSTDKQRLEESQRQFEHILKNEHIKNVPVVL LANKQDMPGALTAEDITRMFKVKKLCSDRNWYVQPCCALTGEGLAQGFRKLTGFVKSHMK SRGDTLAFFKQN*
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	ARL14; ADP-ribosylation factor-like 14; ARF7; ADP-ribosylation factor-like protein 14; ADP-ribosylation factor 7;
-----------------	---

GENE INFORMATION

Entrez Gene ID	80117
-----------------------	-----------------------

UniProt ID	Q8N4G2
-------------------	------------------------

Function	GTP binding; nucleotide binding
-----------------	---------------------------------
