



Mouse anti-Human ARMC4 monoclonal antibody, clone 6G2 (CABT-B9794)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	ARMC4 (NP_060546, 945 a.a. ~ 1045 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	6G2
Conjugate	Unconjugated
Applications	WB,sELISA,ELISA
Sequence Similarities	AISRCCMWGRNRVAFGEHKAVAPLVRYLKSNDTNVHRATAQALYQLSEDADNCITMHENG AVKLLLDVMVGSPDQDLQEAAAGCISNIRRLALATEKARYT*
Format	Liquid
Size	50 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords ARMC4; armadillo repeat containing 4; CILD23; armadillo repeat-containing protein 4;

GENE INFORMATION

Entrez Gene ID	55130
UniProt ID	Q5T2S8
Function	binding