



Mouse anti-Human ARMED monoclonal antibody, clone 3E0 (CABT-B9796)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	ARMED (NP_054873, 287 a.a. ~ 386 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	3E0
Conjugate	Unconjugated
Applications	WB,IP,sELISA,ELISA
Sequence Similarities	AETLAYLIEPDVELQRIASITDHLIAMLADYFKYPSSVSAITDIKRLDHLKHAHEL RQA AFKLYASLGANEDIRKKVSLGEGRPVLTASRQGVST*
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	ARMED; armadillo repeat containing 8; GID5; VID28; S863-2; HSPC056; armadillo repeat-containing protein 8; GID complex subunit 5, VID28 homolog;
-----------------	--

GENE INFORMATION

Entrez Gene ID	25852
UniProt ID	Q8IUR7
Function	binding