



Mouse anti-Human ASZ1 monoclonal antibody, clone 4D0 (CABT-B9806)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	ASZ1 (NP_570124, 127 a.a. ~ 235 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2b
Source/Host	Mouse
Species Reactivity	Human
Clone	4D0
Conjugate	Unconjugated
Applications	WB,sELISA,ELISA
Sequence Similarities	QILKCVELLLSRNADPNVACRRMLTPIMYAARDGHTQVVALLVAHGAEVNTQDENGYTAL TWAARQGHKNIVLKLLELGANKMLQTKDGKMPSEIAKRKNKHHEIFNLL*
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	ASZ1; ankyrin repeat, SAM and basic leucine zipper domain containing 1; ALP1; GASZ; Orf3; ANKL1; C7orf7; CT1.19; ankyrin repeat, SAM and basic leucine zipper domain-containing
-----------------	---

protein 1; ankyrin-like 1; ankyrin-like protein 1; germ cell-specific ankyrin, SAM and basic leucine zipper domain-containing protein;

GENE INFORMATION

Entrez Gene ID	136991
----------------	------------------------

UniProt ID	Q8WWH4
------------	------------------------

Function	signal transducer activity
----------	----------------------------
