



Mouse anti-Human BARHL1 monoclonal antibody, clone 2C8 (CABT-B9836)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	BARHL1 (NP_064448, 55 a.a. ~ 155 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	Human
Clone	2C8
Conjugate	Unconjugated
Applications	WB,sELISA,ELISA
Sequence Similarities	GRDCLETGTPRPGGASGPGLDShLQPGQLSAPAQSRTVTSSFLIRDILADCKPLAACAPY SSSGQPAAPEPGGRLAAKAAEDFRDKLDKSGSNASSDSEY*
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords BARHL1; BarH-like homeobox 1; barH-like 1 homeobox protein;

GENE INFORMATION

Entrez Gene ID	56751
UniProt ID	Q9BZE3
Function	protein binding; sequence-specific DNA binding; sequence-specific DNA binding transcription factor activity