



Mouse anti-Human BCL9L monoclonal antibody, clone 6F9 (CABT-B9846)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	BCL9L (NP_872363, 606 a.a. ~ 716 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Human
Clone	6F9
Conjugate	Unconjugated
Applications	sELISA,ELISA
Sequence Similarities	NQIQRVPGFGGMQSMPEVPMNAMQRPVRPGMGWTEDLPPMGGPSNFAQNTMPYPGGQGE AERFMTPRVREELLRHQLLEKRSMGMQRPLGMAGSGMGQSMEMERMMQAH*
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	BCL9L; B-cell CLL/lymphoma 9-like; BCL9-2; DLNB11; B-cell CLL/lymphoma 9-like protein; protein BCL9-2; BCL9-like protein; B-cell lymphoma 9-like protein; nuclear co-factor of beta-
-----------------	--

catenin signalling;

GENE INFORMATION

Entrez Gene ID	283149
----------------	------------------------

UniProt ID	Q86UU0
------------	------------------------

Function	transcription coactivator activity
----------	------------------------------------
