



Mouse anti-Human BLZF1 monoclonal antibody, clone 3F21 (CABT-B9858)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	BLZF1 (NP_003657, 301 a.a. ~ 401 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	Human
Clone	3F21
Conjugate	Unconjugated
Applications	WB,sELISA,ELISA
Sequence Similarities	KLAKAVNSHLLGNVGINNQQKIPSTVEFCSTPAEKMAETVLRILDPTCKESSPDNPFFE SSPTTLLATKKNIGRFHPYTRYENITFNCCNHCRGELIAL*
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Golgin-45 is a protein that in humans is encoded by the BLZF1 gene. Mouse monoclonal antibody raised against a partial recombinant BLZF1.
---------------------	---

Keywords BLZF1; basic leucine zipper nuclear factor 1; JEM1; JEM-1; JEM-1s; GOLGIN-45; Golgin-45; JEM-1short protein; cytoplasmic protein; p45 basic leucine-zipper nuclear factor;

GENE INFORMATION

Entrez Gene ID [8548](#)

UniProt ID [Q9H2G9](#)

Function DNA binding; enzyme binding; sequence-specific DNA binding transcription factor activity; ubiquitin protein ligase binding
