



Mouse anti-Human BUD31 monoclonal antibody, clone 2F21 (CABT-B9870)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	G10 (NP_003901, 35 a.a. ~ 145 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2b
Source/Host	Mouse
Species Reactivity	Human
Clone	2F21
Conjugate	Unconjugated
Applications	WB,sELISA,ELISA
Sequence Similarities	EPHEGKRKVESLWPIFRIHHQKTRYIFDLFYKRKAISRELYEYCIKEGYADKNLIAKWKK QGYENLCCLRCIQTRDTNFGTNCICRVPKSKLEVGRRIIECTHCGCRGCSG*
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	BUD31; BUD31 homolog (S. cerevisiae); G10; EDG2; Cwc14; EDG-2; fSAP17; YCR063W; protein BUD31 homolog; protein EDG-2; protein G10 homolog; maternal G10 transcript; G10
-----------------	---

GENE INFORMATION

Entrez Gene ID	8896
UniProt ID	P41223
Pathway	Spliceosome, organism-specific biosystem; Spliceosome, conserved biosystem
Function	sequence-specific DNA binding transcription factor activity