



Mouse anti-Human C20orf7 monoclonal antibody, clone 5B0 (CABT-B9873)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	C20orf7 (AAH05984, 1 a.a. ~ 159 a.a) full length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2b
Source/Host	Mouse
Species Reactivity	Human
Clone	5B0
Conjugate	Unconjugated
Applications	ELISA
Sequence Similarities	MFGGDTLYELRCSLQLAETEREGGFSPHISPFTAVNDLGHLLGRAGFNTLTVDTDEIQVN YPGMFELMEDLQGMGESNCAWNRKALLHRDTMLAAAAYREMYRNEDGSVPATYQIYYMI GWKYHESQARPAERGSATVSFGELGKINNLMPPGKKSQ*
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords chromosome 20 open reading frame 7; bA526K24.2; dJ842G6.1; FLJ22324; EC 2.1.1.-;

MGC90272; EC 2.1.1; probable methyltransferase C20orf7, mitochondrial;
OTTHUMP00000030304; OTTHUMP00000030305; OTTHUMP00000162656

GENE INFORMATION

Entrez Gene ID	79133
----------------	-----------------------

UniProt ID	Q5TEU4
------------	------------------------

Function	methyltransferase activity; molecular_function; transferase activity
----------	--
