



Mouse anti-Human C9orf96 monoclonal antibody, clone 3C6 (CABT-B9876)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	C9orf96 (NP_714921, 1 a.a. ~ 100 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	Human
Clone	3C6
Conjugate	Unconjugated
Applications	WB,sELISA,ELISA,RNAi Knockdown
Sequence Similarities	MLGPGSNRRRPTQGERGPGSPGEPMEKYQVLYQLNPGALGVNLVVEEMETKVKHVIKQVE CMDDHYASQALEELMPLLKLRHAHISVYQELFITWNGEI*
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords chromosome 9 open reading frame 96; SgK071; MGC43306; Sk521; Sugen kinase 071; protein kinase-like protein SgK071; RP11-244N20.8; OTTHUMP00000022465

GENE INFORMATION

Entrez Gene ID	169436
UniProt ID	Q8NE28
Function	ATP binding; nucleotide binding; protein kinase activity; transferase activity, transferring phosphorus-containing groups