



Mouse anti-Human CARD16 monoclonal antibody, clone 4F21 (CABT-B9900)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	COPI (NP_443121, 1 a.a. ~ 98 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2b
Source/Host	Mouse
Species Reactivity	Human
Clone	4F21
Conjugate	Unconjugated
Applications	WB,sELISA,ELISA
Sequence Similarities	MADKVLKEKRKLFHSMGEGTINGLLDELLQTRVLNQEEMEKVKRENATVMDKTRALIDS VIPKGAQACQICITYICEEDSYLAETLGLSAGPIPGN*
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	CARD16; caspase recruitment domain family, member 16; COP; COP1; PSEUDO-ICE; caspase recruitment domain-containing protein 16; CARD only protein; CARD-only protein 1;
-----------------	---

caspase-1 inhibitor COP; pseudo-IL1B-converting enzyme; CARD only domain-containing protein 1; caspase recruitment domain-only protein 1; pseudo interleukin-1beta converting enzyme; pseudo interleukin-1 beta converting enzyme; caspase-1 dominant-negative inhibitor pseudo-ICE;

GENE INFORMATION

Entrez Gene ID	114769
----------------	------------------------

UniProt ID	Q5EG05
------------	------------------------

Function	cysteine-type endopeptidase inhibitor activity; peptidase inhibitor activity
----------	--
