



Mouse anti-Human CCDC132 monoclonal antibody, clone 3E22 (CABT-B9906)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	FLJ20097 (NP_060137, 862 a.a. ~ 964 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 37KDa.
Isotype	IgG2b
Source/Host	Mouse
Species Reactivity	Human
Clone	3E22
Conjugate	Unconjugated
Applications	WB, IHC, IF
Sequence Similarities	GYANVKKCSNEGRALMQLDFQQFLMKLEKLTDIRPIPDKEFVETYIKAYYL TENDMERWI KEHREYSTKQLTNLVNVCLGSHINKKARQKLLAAIDDIRPKR
Format	Liquid
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords CCDC132; coiled-coil domain containing 132; coiled-coil domain-containing protein 132;

GENE INFORMATION

Entrez Gene ID [55610](#)

UniProt ID [Q96JG6](#)
