



Mouse anti-Human CCNL1 monoclonal antibody, clone 2G82D6 (CABT-B9912)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	CCNL1 (AAH38394, 1 a.a. ~ 173 a.a) full length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	2G82D6
Conjugate	Unconjugated
Applications	WB,IF,sELISA,ELISA
Sequence Similarities	MASGPHSTATAAAAASSAAPSAGGSSSGTTTTTTTTTGGILIGDRLYSEVSLTIDHSLIP EERLSPTPSMQDGLDLPSETDLRILGCELIQAAGILLRLPQVAMATGQVLFHRFFYSKSF VKHSFEIVAMACINLASKIEEAPRRIRDVINVFHHLRQLRGKSDQLHLPKPG*
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction Cyclin-L1 is a protein that in humans is encoded by the CCNL1 gene. Mouse monoclonal

antibody raised against a full length recombinant CCNL1.

Keywords

CCNL1; cyclin L1; ANIA6A; BM-001; PRO1073; ania-6a; cyclin-L1; cyclin-L; cyclin L gamma; cyclin L ania-6a;

GENE INFORMATION

Entrez Gene ID

[57018](#)

UniProt ID

[Q9UK58](#)

Function

protein kinase binding
