



Mouse anti-Human CCT8L2 monoclonal antibody, clone 5H7 (CABT-B9914)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	CESK1 (AAH33797, 1 a.a. ~ 558 a.a) full length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	Human
Clone	5H7
Conjugate	Unconjugated
Applications	WB,ELISA,RNAi Knockdown
Sequence Similarities	MDSTVPSALELPQRLALNPRESPRSPEEEEPHLLSSLAAVQTLASVIRPCYGPBGRQKFL VTMKGETVCTGCATAILRALELEHPAAWLLREAGQTQAENSGDGTAFVVLLTEALLEQAE QLLKAGLPRPQLREAYATATAEVLATLPSLAIQSLGPLEDPSWALHSVNMNTHLSPMDHL TKLVAHACWAIKELDGSFKPERVGVICALPGGTLEDSCLLPGLAISGKLCGQMATVLSGAR VALFACPFPGPAHPNA
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	CCT8L2; chaperonin containing TCP1, subunit 8 (theta)-like 2; CESK1; putative T-complex protein 1 subunit theta-like 2; T-complex protein 1;
-----------------	--

GENE INFORMATION

Entrez Gene ID	150160
-----------------------	------------------------

UniProt ID	Q96SF2
-------------------	------------------------

Function	ATP binding; anion channel activity; calcium-activated potassium channel activity; nucleotide binding
-----------------	---
