



Mouse anti-Human CDC20B monoclonal antibody, clone 3G3 (CABT-B9919)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	CDC20B (NP_689836, 416 a.a. ~ 515 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	3G3
Conjugate	Unconjugated
Applications	WB,IF,sELISA,ELISA
Sequence Similarities	VLAIGGGMKDGRHLHILDINAGKSIQTPSTNSQICSLIWLPKTKEIATGQGTPKNDVTVWT CPTVSRSGHRGRVLHLALSPDQTRVFSAAADGTASVWNCY
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	CDC20B; cell division cycle 20B; G6VTS76519; cell division cycle protein 20 homolog B; CDC20-like protein; cell division cycle 20 homolog B; CDC20 cell division cycle 20 homolog B;
-----------------	---

GENE INFORMATION

Entrez Gene ID [166979](#)

UniProt ID [Q86Y33](#)
