



Mouse anti-Human CDKN2AIP monoclonal antibody, clone 2B7 (CABT-B9933)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	CARF (AAH22270,2a.a. ~ 579 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	2B7
Conjugate	Unconjugated
Applications	WB, IHC, IF, sELISA, ELISA
Sequence Similarities	MAQEVSEYLSQNPRVAAWVEALRCDGETDKHWRHRRDFLLRNAGDLAPAGGAASASTDEA ADAESGTRNRQLQQILISFSMAWANHVFLGCRYPQKVMDKILSMAEGIKVTDAPTYTTRDE LVAKVKKRGISSNEGVEEPSKKRVIEGKNSSAVEQDHAKTSAKTERASAQQENSSTCIG SAIKSESGNSARSSGISSQNSSTSDGDRSVSSQSSSSVSSQVTTAGSGKASEAEAPDKHG SSFVSLKSSVN
Format	Liquid
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords CDKN2AIP; CDKN2A interacting protein; CARF; CDKN2A-interacting protein; collaborator of ARF; collaborates/cooperates with ARF (alternate reading frame) protein;

GENE INFORMATION

Entrez Gene ID [55602](#)

UniProt ID [Q9NXV6](#)

Function double-stranded RNA binding; p53 binding; protein binding;
