



Mouse anti-Human CFDP1 monoclonal antibody, clone 6C8 (CABT-B9954)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	CFDP1 (NP_006315, 168 a.a. ~ 252 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	6C8
Conjugate	Unconjugated
Applications	WB,ELISA
Sequence Similarities	VFDFAGEEVRVTKEVDATSKEAKSFFKQNEKEKPQANVPSALPSLPAGSGLKRSSGMSSL LGKIGAKKQKMSTLEKSKLDWESF*
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Craniofacial development protein 1 is a protein that in humans is encoded by the CFDP1 gene. Mouse monoclonal antibody raised against a partial recombinant CFDP1.
---------------------	--

Keywords

CFDP1; craniofacial development protein 1; p97; BCNT; CP27; SWC5; Yeti; CENP-29; BUCENTAUR; centromere protein 29; phosphoprotein (Bucentaur);

GENE INFORMATION

Entrez Gene ID	10428
UniProt ID	Q9JEE9
Function	molecular_function