



Mouse anti-Human CGRRF1 monoclonal antibody, clone 3F23 (CABT-B9958)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	CGRRF1 (AAH15063, 36 a.a. ~ 333 a.a) full length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	3F23
Conjugate	Unconjugated
Applications	sELISA,ELISA
Sequence Similarities	MAAVFLVTLYEYSPLFYIAVVFTCFIVTTGLVLGWFGWDVPVILRNSEETQFSTRVFKKQ MRQVKNPFGLEITNPSSASITTGITLTDCLEDSLLTCYWGCNVQKLYEALQKHVYCFRI STPQALEDALYSELYQEYFIKKDSKEEICYQLPRDTKIEDFGTVPRSRYPVALLTLA DEDDREIYDIISMVSVIHIPDRTYKLSCRILYQYLLAQGFHDLKQLFMSANNFTPSN NSSSEEKNTDRSLLE
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	CGRRF1; cell growth regulator with ring finger domain 1; CGR19; RNF197; cell growth regulator with RING finger domain protein 1; RING finger protein 197; cell growth regulatory gene 19 protein;
-----------------	---

GENE INFORMATION

Entrez Gene ID	10668
-----------------------	-----------------------

UniProt ID	Q99675
-------------------	------------------------

Function	metal ion binding; zinc ion binding
-----------------	-------------------------------------
