



Mouse anti-Human CHERP monoclonal antibody, clone 3I6 (CABT-B9962)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	CHERP (NP_006378, 795 a.a. ~ 884 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2b
Source/Host	Mouse
Species Reactivity	Human
Clone	3I6
Conjugate	Unconjugated
Applications	WB, ELISA
Sequence Similarities	GSNSAPPIPDSRLGEENKGHQLVKMGWSGSGGLGAKEQGIQDPIKGGDVRDKWDQYKGV GVALDDPYENYRRNKSYSFIARMKARDEC*
Format	Liquid
Size	50 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	CHERP; calcium homeostasis endoplasmic reticulum protein; SRA1; DAN16; SCAF6; ERROT213-21; ERROT 213-21; protein with polyglutamine repeat; SR-related CTD
-----------------	--

associated factor 6; SR-related CTD-associated factor 6;

GENE INFORMATION

Entrez Gene ID	10523
----------------	-----------------------

UniProt ID	Q8IWX8
------------	------------------------

Pathway	Spliceosome, organism-specific biosystem; Spliceosome, conserved biosystem
---------	--

Function	RNA binding; molecular_function
----------	---------------------------------
