



Mouse anti-Human CHML monoclonal antibody, clone 6H5 (CABT-B9963)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	CHML (NP_001812, 68 a.a. ~ 178 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	Human
Clone	6H5
Conjugate	Unconjugated
Applications	WB,sELISA,ELISA
Sequence Similarities	IGESTVWQDLIHETEEAITLRKKDETIQHTEAFQYASQDMEDNVEEIGALQKNPSLGV SNTFTEVLDSALPEESQLSYFNSDEMPAKHTQKSDTEISLEVTDVEESVE*
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	The product of the CHML gene supports geranylgeranylation of most Rab proteins and may substitute for REP-1 in tissues other than retina. CHML is localized close to the gene for Usher
---------------------	---

syndrome type II. [provided by RefSeq, Jul 2008]

Keywords

CHML; choroideremia-like (Rab escort protein 2); REP2; rab proteins
geranylgeranyltransferase component A 2; choroideremia-like protein; choroideraemia-like
protein;

GENE INFORMATION

Entrez Gene ID

[1122](#)

UniProt ID

[P26374](#)

Function

GTPase activator activity; Rab geranylgeranyltransferase activity
