



Mouse anti-Human CNIH3 monoclonal antibody, clone 2F21 (CABT-B9994)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	CNIH3 (AAH22780, 1 a.a. ~ 161 a.a) full length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2b
Source/Host	Mouse
Species Reactivity	Human
Clone	2F21
Conjugate	Unconjugated
Applications	ELISA
Sequence Similarities	MAFTFAAFCYMLSLVLCAALIFFAIWHIIAFDELRTDFKSPIDQCNPVHARERLRNIERI CFLLRKLVLP EYSIHSLFCIMFLCAQEWLTGLNVPLL FYHFWRYFHCPADSSELAYDPP VVMNADTLSYCQKEAWCKLAFYLLSFFYYLYCMIYTLVSS*
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords CNIH3; cornichon homolog 3 (Drosophila); CNIH-3; 2900075G08Rik; protein cornichon

GENE INFORMATION

Entrez Gene ID [149111](#)

UniProt ID [Q8TBE1](#)
