



Mouse anti-Human CNIH4 monoclonal antibody, clone 5F212H5 (CABT-B9995)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	HSPC163 (AAH00573, 1 a.a. ~ 140 a.a) full length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	Human
Clone	5F212H5
Conjugate	Unconjugated
Applications	WB,ELISA
Sequence Similarities	MEAVVFVFSLLDCCALIFLSVYFIITLSDLECDYINARSCCSKLNKWWIPELIGHTIVTV LLLMSLHWFI LLNLPVATWNIYRYIMVPSGNMGVFDPTTEIHNRGQLKSHMKEAMIKLGF HLLCFFMYLYSMILALIND*
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords CNIH4; cornichon family AMPA receptor auxiliary protein 4; CNIH-4; HSPC163; protein

GENE INFORMATION

Entrez Gene ID	29097
----------------	-----------------------

UniProt ID	Q9P003
------------	------------------------

Function	protein binding
----------	-----------------
